

FOX

Diode laser



PORTABLE - PRACTICAL - POWERFUL

- ▶ Integrated Power Meter for Calibration
- ▶ Video Player Feature
- ▶ Green Aiming Beam
- ▶ Touch Screen Display
- ▶ Low Level Therapy Mode
- ▶ Up to 11 W Distal Output Power



A.R.C.
LASER


enlighten your surgery.



calibrate the FOX using the internal POWER METER.




Watch the latest training movie on the screen.



the green aiming beam is highly visible.



Modern touch Screen Display



portable and mobile - ideal for outpatient centres.

FOX CALIBRATION

FOX is a leader in performance.

- Connect directly to the FOX sensor via USB to measure, calibrate and verify laser power output at the fiber distal tip.

FOX VIDEO

FOX is a leader in innovation.

- Access and view the A.R.C. library of videos for demonstration and training of the latest FOX laser applications and techniques.

GREEN AIMING BEAM

FOX is a leader in design.

- NEW green aiming beam enhances visibility and enables precise positioning of the laser for treatment.

TOUCH SCREEN DISPLAY

FOX is a leader in advanced features.

- Intuitive Touch Screen Display puts quick and easy treatment parameter selections at your fingertip.

FOX PORTABILITY

FOX is a leader in freedom of use.

- Battery powered operation enables FOX laser usage anytime, anywhere using the convenient carrying case.

ENT

Ambulatory Treatment: Endonasal Surgery

- Mucosal Hyperplasia
- Polyposis
- Cysts, mucous cysts
- Epistaxis, Morbus Osier
- Hemangiomas, Telangectasls
- Synechia, Stenosis in endonasal Structures
- Septal Spur, Septum deformation
- Nasal sinus surgery

Dacryocystorhinostomy (DCR)

Oro-pharynx

- Uvulopatectomy (LAUP)
- Tonsillotomy, Tonsillectomy
- Tumor Vaporization

Stationary Treatment:

- Otology
- Stapedectomy / Stapedotomy
- Removal of Cholesteatoma
- Myringotomy

Larynx

- Vocal Polypctomy
- Papilloma
- Laryngeal Carcinoma
- Pediatrics
- Vocal nodules
- Congenital Laryngeal Stenosis
- Neoplastic Stenosis
- Tumovaporisation Laser assisted
- Tuboplasty of the eustachian Tube (LETP) with an special Micro-endoscopy

Endoscopy / Microendoscopy

- Eustachian Tube
- Lacrimal Passages
- Salivary Endoscopy
- Micro-Otology
- Micro Sinus Surgery

DENTISTRY

- Endodontics
- Germ reduction in the root channel, also inside the dentintubuli
- Parodontology
- Gingival Treatment
- Periimplantitis
- Gingival forming
- Oeantng of the Gingival Bags and Germ reduction
- Dental Surgery
- Lip-/lingual mucous membrane
- Covered Implants and tooth
- Coagulation
- Herpes
- Cysts
- Fistula walks
- Fibrosis
- Prosthesis mould areas

DERMATOLOGY/ AESTHETICS

- Sclerosis
- Coagulation
- Skyrocket capillary ectasis, spider-burst
- Vascular Lesions
- Teteangioma
- Spider Veins
- Hemangioma

SURGERY

- Coagulation
- Interstitial Surgery
- Minimal invasive Intrusions
- Treatment of vascular lesions surface vessels
- Sderotherapy of widened veins




Work sterile with our protective cap.

FOX SAFETY

FOX is a leader in safety.

- Sterile, protective caps for use on the "I-drive" control knob creates safe and easy use in a sterile field.

FOX SETS NEW LASER STANDARDS OF THE FUTURE -TODAY !



FOX is the most PORTABLE, PRACTICAL & POWERFUL laser of its kind and offers unmatched VALUE for a broad range of laser applications in multiple treatment settings. FOX is leading the way.

Specifications

Parameter	Diode
Wavelength, output power*	810 nm 8 Watt 940 nm 4 Watt 980 nm 11 Watt 1064 nm 10 Watt
Pulse length	1 ms to cw
Pulse intervall	4 ms to 1 sec.
Aiming beam	Green 532 nm, < 2 mW
Dimensions (W x D x H)	12 x 21 x 10 cm
Weight	1.2 kg ready to use
Power requirement	Rechargeable batteries

wavelength at order placement. Specifications are subject to change without notice



Marketed By :

Medilife Technologies
1st Floor, Plot No.2, Vikash Plaza Local Shopping
Above HDFC Bank Kalkaji New Delhi-110019
Ph: +91 11 40588606/07, Mob: 9810524302
Telefax : 011-40588606
Email: info@medilife.in, sales@medilife.in

Manufactured By :

A.R.C.
LASER

